




ENVIRO



CENTRIFUGE TUBES FOR PAH EXTRACTION



ENVIRO-CLEAN® PAH

UCT's ENVIRO-CLEAN® PAH certified centrifuge tubes ECPAHFR50CT (C of A available) are specially designed for performing PAH analysis using either QuEChERS, AOAC or other methods that require the use of 50 mL centrifuge tubes. UCT's special polypropylene tubes with plug-seal caps are ideal for performing low level PAH extractions.

	Specification Values	Analysis
Acenaphthene	0.5 ppb	ND
Acenaphthylene	1 ppb	< 0.5 ppb
Anthracene	0.5 ppb	ND
Benz(a)anthracene	0.5 ppb	ND
Benzo(a)pyrene	0.5 ppb	ND
Benzo(b)fluoranthene	0.5 ppb	ND
Benzo(g,h,i)perylene	0.5 ppb	ND
Benzo(k)fluoranthene	0.5 ppb	ND
Chrysene	0.5 ppb	ND
Dibenzo(a,h)anthracene	0.5 ppb	ND
Fluoranthene	0.5 ppb	ND
Fluorene	0.5 ppb	ND
Indeno(1,2,3-cd)pyrene	0.5 ppb	ND
Naphthalene	1 ppb	< 0.3 ppb
Phenanthrene	0.5 ppb	ND
Pyrene	0.5 ppb	ND

*Note : Values based on 10g sample size

Placing An Order

Phone: 215.781.9255
800.385.3153
Fax: 215.785.1226

UCT, Inc.
2731 Bartram Rd.
Bristol, PA 19007

Email: info@unitedchem.com
Web: www.unitedchem.com



4101-03-01

